

Treatment of Asian Skin with an Advanced Intense Pulsed Light Source

Hartne Suzuki, MD

Suzuki Plastic Surgery Clinic, Kyoto Japan

Summary

Intense Pulsed Light (IPL) sources have been improved to the point that they rival multicolor pulsed lasers in the visible and near-infrared spectrum for treatment of many vascular and pigmentation disorders. Compared with second generation IPLs, earlier IPLs had a broader spectral band including infrared wavelengths absorbed by water and other smaller pigment targets. Asian skin is particularly prone to post-inflammatory hyperpigmentation, and post-laser IPL (Vascular IPL™, ISC-Lang PL™, 580 pulses, 300 ms pulse width, 1.3-2.8 J/cm²; 2.4 ms and 5 ms pulses with 30 ms delay) for a second generation IPL (Ellipse-Flux™, ISC-Lang PL™, 1000-700 pulses; V1.2 or P1 handpieces; 9.12 J/cm² or 7.9 J/cm²; 2.5 ms and 2.5 ms pulses with 30 ms delay). Most of the patients (66%) were treated for facial lesions associated with photodamage including lentigo, melasma, and fine wrinkles. Other patients (94%) were treated for congenital nevi, scars, or periorbital lesions. Efficacy was comparable with both IPLs, but was obtained with much lower fluence (1.3), fewer treatments (1.2) and fewer side effects (1.7) using the second generation IPL. Side effects occurred in 7% of patients after the first generation IPL and 1% with the second generation IPL, consisting of blistering, inflammation and transient pigmentary changes. There was no scarring. In summary, there are strong advantages of second generation IPL.

Purpose

The purpose for the study was to compare the efficacy and side effects of the first and the second generation IPL systems.

Equipments and Parameters

The first IPL system, Vascular IPL™, ISC-Lumant, 580 nm or 590 nm filter handpieces; 31-35 J/cm²; 2.4 ms and 3 ms pulse with 30 ms delay.

The second IPL system, Ellipse-Flux™, ISC, V1.2 or P1 handpieces; 9-12 J/cm² or 7-9 J/cm²; 2.5 ms and 2.5 ms pulse with 30 ms delay.

The main difference between the first and the second IPL system is that the second generation IPL system has dual mode filtering, the single mode filtering cut-off only shorter wavelengths including ultraviolet, whereas the dual mode filtering also eliminates longer wavelengths which are absorbed mainly by water, which unnecessarily heats the skin. Dual mode filtering enables most of the energy applied to the skin to work effectively for hemoglobin or melanin treatment.

Methods

For rejuvenation entire face applying the upper eyelids which were covered with metal eye shields, were treated with the parameters shown above.

Patient with pigmentation problems were treated mainly by 580 nm filter handpieces with the Vascular IPL™ and V1.2 handpiece with the Ellipse-Flux™.

Patient with vascular problems were mainly treated with 590 nm filter handpieces with the Vascular IPL™ and P1 handpiece with the Ellipse-Flux™.

Treatment was repeated every 3 weeks until the patients were satisfied or no further improvements were seen.

Patients were allowed to make up next day of the treatments.

Results


Similar endpoint was obtained with 31-35 J/cm² on with the Vascular IPL™ and 7-12 J/cm² on with the Ellipse-Flux™ and with less pain.

Five hundred cases including 489 rejuvenation cases were treated with the Vascular IPL™ from April 2000 to June 2002. Six hundred and 99 cases including 648 rejuvenation cases were treated with the Ellipse-Flux™ from July 2002 to March 2003. The average number of treatments for rejuvenation were 7.2 for the Vascular IPL™ and 4.2 for the Ellipse-Flux™. Side effects were seen in 7% of the cases treated with the Vascular IPL™. Side effects with the Ellipse-Flux™ were seen only when we had repeatedly using the IPL or IPL applicator for restricted areas. Ellipse-Flux™ for rejuvenation using typical parameter never have caused any side effects.

Conclusion

A second generation IPL system, Ellipse-Flux™ using the dual mode filtering, requires fewer treatments with lower fluence for rejuvenation on the face, compared with the first generation IPL. The second generation IPL mainly causes side effects, or downtime.

Difference between Single and Dual Mode Filtering



1) The diagram shows the wavelength emitted from the first generation IPL system. You have to use wavelengths in this range to treat pigmentation problems. However, shorter wavelengths are readily absorbed by water.


2) The second generation IPL system has a dual mode filtering, which eliminates both shorter and longer wavelengths. The shorter wavelengths are filtered out by the IPL applicator we use for most of the cases for rejuvenation.

3) The wavelength through the IPL applicator, which is 580 to 790 nm, covers the spectral range for hemoglobin and melanin. Shorter wavelengths are not absorbed by vascular lesions.

Conditions treated with each IPL




Side effects




Rejuvenation cases with the second generation IPL system


Facial pigmented lesion




Elevated pigmented lesion



Melasma and solar lentiginies



Moles



Facial erythema and wrinkles

