



SIRIUS



The Fastest and Most Cost-Effective Hair Removal and Skin Rejuvenation System

HIGH SPEED

With its high speed and superior technology, Sirius reduces treatment time while at the same time increasing the efficacy of your treatments.

SAFE AND EFFECTIVE TREATMENTS

With Sirius you simply select the treatment, the patient's skin type, suntan, vessel size or hair thickness. Sirius then automatically chooses the appropriate settings including the energy level and pulse length. Ellipse's clinical intelligence ensures safe and effective treatments and reduces the risk of side effects. Moreover, a built-in patient database stores all client and treatment data for you.

RANGE OF TREATMENTS

Sirius comes with indications for use on hair removal, facial thread veins and skin rejuvenation.

FAST AND RELIABLE

With the all new Remote Assistance feature you get the markets most reliable and state-of-the-art online system. Thanks to a secure online access your system can be diagnosed, updated and even get the latest software and newest applications. This can be done whenever it suits you to connect - even after work hours.

Safe and effective treatments

Fast treatments and high client throughput

Efficient and comfortable SWT[®] applicators



Silent and ergonomic operation

Remote Assistance




EXCELLENT OPERATOR COMFORT

With the SWT[®] applicators you reduce the risk of repetitive strain injury. The light weight of the applicators and the soft and flexible cables make it easy to change position. The design of the applicator head gives you a clear view of the treatment area in a comfortable working position, and for your convenience the accessory tray with screen can be turned to the left and right. Moreover, Sirius is the most silent system on the market, and both you and the patient can enjoy a calm and comfortable working environment.




SWT® APPLICATORS

You can connect two lightweight Selective Waveband Technology (SWT®) applicators at the same time and even change applicators while the system is running. For Sirius, Ellipse's exclusive SWT® applicators come in 5 narrowband versions. Among the versions included are the HRD 645 for hair removal on dark skin types 4-6 and the large HRL 600 applicator with an astonishing 8.64 cm² footprint. The large footprint of the HRL 600 allows you to treat a full back in less than 10 minutes, and fewer overlaps and fewer skipped areas ensure a much better result. For rejuvenation treatments, you also have the choice between the VL 555 and the VLS 555 – a small spot applicator which allows you to treat close to and around the eyes, and around the wings of the nose, in corners and folds.



EFFECTIVE TREATMENTS WITH HIGH PATIENT COMFORT


With SWT® technology, unnecessary and potentially harmful wavelengths are eliminated, allowing only those that are beneficial to individual treatments to reach the skin. This significantly reduces the amount of energy required, so improving patient comfort. Combined with perfect square pulse technology, and Ellipse patented dual mode filter, SWT® enables you to treat skin types 1-6 without the need for active skin cooling. The result is safer, more predictable and more effective treatments.



Before & After



Skin Rejuvenation
1 Tx (VL 555)
Courtesy of Ellipse Spain



Hair Removal (skin type 5)
5 Tx (HR 600)
Courtesy of Dr. Khan



Facial Thread Veins
3 Tx (VL 555)
Courtesy of Dr. Sandra Smaliene



Mieke Morgan, Laser Skin Specialist – The London Cosmetic Clinic, London, United Kingdom

“Both the VL 555 and the HR applicators have been designed with the operator in mind. The quick charging time in between shots, the swivel base, soft flexible cord and defined crystal block allow for easy access to all areas of the body as well as fast, successful treatments.



The HR comes in two sizes allowing treatment times to be reduced drastically. I enjoy using the standard size for smaller areas such as the face, underarms etc. And the larger applicator for the body. The specific settings associated with these applicators have improved my results considerably and I find treating stubborn patches or finer hairs much easier.

The VL 555 is possibly my favourite applicator. It truly is versatile and allows for many concerns to be addressed as well as many skin treatments to be performed. My clients have commented on how much more comfortable the new applicators are as well as how much less the stingy, red downtime is. I found it to be a great investment that really has expanded my treatment portfolio.”

FLEXIBLE PARTNERSHIP

We would like to welcome you into the world of Ellipse – where perfection, technical superiority, quality and long-term partnership meet. With Sirius we offer you low running costs and no hidden ownership cost. A one-year license lets you fire limitless shots in the first year at no extra expense.




Sirius is a trademark of Ellipse A/S, a producer of light-based solutions for medical and cosmetic conditions. Ellipse A/S has gained the acclaim of doctors worldwide for the quality of Ellipse products and services.

Ellipse A/S. All rights reserved.



Ellipse A/S • Agern Allé 11 • DK-2970 Hørsholm, Denmark • ellipse.com



ISO 9001 ISO 13485 CE 0543

