

Nordlys: Indications for Use

	Skin type						Applicators
	1	2	3	4	5	6	
SWT TREATMENTS							
Hair Removal	●	●	●	●	●		HR 600 / HRL 600
Hair Removal	●	●	●	●	●	●	HRD 645
Photorejuvenation (Diffuse Redness and Pigment)	●	●	●	●			PR(S) 530 / VL(S) 555
Telangiectasias Face	●	●	●	(●)			PR(S) 530 / VL(S) 555
Port Wine Stains	●	●	●	(●)			PR(S) 530 / VL(S) 555
Acne	●	●	●	●			PR(S) 530
Rosacea	●	●	●	(●)			PR(S) 530 / VL(S) 555
Poikiloderma of Civatte	●	●	●	(●)			PR(S) 530 / VL(S) 555
Skin Texture*	●	●	●	●			PR(S) 530
PDT Photorejuvenation*	●	●	●	●	●		PL 400 / PR(S) 530
PDT Acne*	●	●	●	●	●		PL 400 / PR(S) 530
Pigmented Lesions	●	●	●	●	●		PL 400

*Treatment of Skin Texture and PDT (Photodynamic Therapy) Assisted Treatments of Acne and Photorejuvenation are not FDA cleared, and are not available in US.

(●) Only VL 555

LASER TREATMENTS							
Leg Veins	●	●	●	●	●	●	Nd:YAG 1064
Telangiectasias and Reticular Vessels	●	●	●	●			Nd:YAG 1064
Venous Lakes	●	●	●	●			Nd:YAG 1064
Port Wine Stains	●	●	●	●			Nd:YAG 1064
Onychomycosis (Clear Nail)	●	●	●	●	●		Nd:YAG 1064
Warts	●	●	●	●	●		Nd:YAG 1064
Acne Scars	●	●	●	●	●		Frax 1550
Surgical Scars	●	●	●	●	●		Frax 1550
Non-ablative Skin Resurfacing	●	●	●	●	●		Frax 1550
Stretch Marks	●	●	●	●	●		Frax 1550



www.ellipse.com

nordlys

MORE THAN LIGHT
by Ellipse

SWT[®] Selective Waveband Technology



Nd:YAG 1064 Laser



Frax 1550 Fractional Non-Ablative Laser

